


EG00700

Two Channel Temperature Module



°C

oem.medlab.eu

- Two independent channels plus reference channel
- Small size (77 x 44 mm), maximum height 10 mm
- Single power supply (5 V DC), 50 mW
- RS232 and TTL level serial asynchronous interface
- Integrated patient isolation (4000 V RMS, IEC 601 CF type)
- Integrated sensor-off detection
- Transmission of temperature values, resolution 0.1 °C
- Measuring range 20.0 °C up to 44.0 °C
- Compatible with all YSI 400 Probes (2252 Ohms at 25 °C)
- Accuracy: +/- 0.1°C + probe tolerance
- Made in Germany

Applications:

Patient Monitors

Temperature measurement for incubators aso.



medlab

medizinische Diagnosegeräte GmbH

Helmholtzstrasse 1a

D 76297 Stutensee near Karlsruhe

Germany

Tel.: +49 (0) 7244 74110 0

FAX: +49(0) 7244 74110 28

sales@medlab.eu