


# EG00700

## Two Channel Temperature Module



[oem.medlab-gmbh.de](http://oem.medlab-gmbh.de)

- Two independent channels plus reference channel
- Small size ( 77 x 44 mm ), maximum height 10 mm
- Single power supply ( 5 V DC ), 50 mW
- RS232 and TTL level serial asynchronous interface
- Integrated patient isolation according to IEC601-1
- Integrated sensor-off detection
- Transmission of temperature values, resolution 0.1 °C
- Measuring range 20.1 °C up to 45.0 °C
- Compatible with all YSI 400 Probes
- Accuracy: +/- 0.1°C + Probe Tolerance
- Made in Germany

Applications:

Patient Monitors

Temperature measurement for incubators aso.



**medlab**

medizinische Diagnosegeräte GmbH

Helmholtzstrasse 1

D 76297 Stutensee near Karlsruhe

Germany

Tel.: +49 (0) 7244 74110 0

FAX: +49(0) 7244 74110 28

[sales@medlab-gmbh.de](mailto:sales@medlab-gmbh.de)