


SpO₂-Probes

High quality probes for
Medlab Pulse Oximeters



Medlab SpO₂ probes can be used exclusively with Medlab pulse oximeters and Medlab SpO₂-Modules. They are the result of more than 15 years of experience in developing SpO₂ sensors.

By using only the highest quality cables, Kevlar reinforcement and cable strain relief, we are able to offer all our reusable probes with a two-year warranty.

We deliver probes for different applications and patient groups.

Sensor extensions are optionally available.



Extension cable DSUB-DSUB
120 cm
Model No. 00310



Extension cable Lemo-DSUB
150 cm
Model No. 00309

Medlab medizinische Diagnosegeräte GmbH
Helmholtzstr. 1a, 76297 Stutensee - Germany
Tel. +49 (0) 7244 74110-0, FAX +49 (0) 7244 74110-28
sales@medlab.eu, www.medlab.eu

medlab
medizinische Diagnosegeräte GmbH

Ordering Info

Reusable Probes

Type	Patient Weight	Model No.
P-200 Fingerclip Probe	> 20 kg	00535
PX-200 Fingerclip Probe	> 20 kg	00536
PS-200 Fingerclip Probe small	10 - 50 kg	00537
R-200 Silicone Probe	> 20 kg	00550
Y-200 Multisite Probe	> 1 kg	00526



P-200
Adult Fingerclip Probe,
140 cm
standard delivery



PX-200
Adult Fingerclip Probe,
280 cm



PS-200
Pediatric Fingerclip Probe,
140 cm



R-200
Adult Oxiflex Silicone Probe,
140 cm



Y-200
Multisite Probe,
120 cm

Disposable Probes - Adhesive Plaster

Type	Applicable for	Model No.
DIS-200AP*	Adult > 20 kg	00540
DIS-200CP*	Child 8 - 30 kg	00541
DIS-200IP**	Infant 5 - 15 kg	00542
DIS-200NP**	Neonate < 5 kg	00543

* Length 45cm

** Length 90cm



Disposable Probes
Delivery:
Package of 24 pieces

CE 0494